**TERRO**

VÜCUT ANALİZ BASKÜLÜ - KULLANIM KILAVUZU  
BODY ANALYSIS SCALE - INSTRUCTION MANUAL

## Giriş

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INTRODUCTION

Thank you for purchasing Fakir Terro Body Analysis Scale which has been designed and manufactured to give you many years of trouble-free service. You may already be familiar with using a similar product, but do please take time to read these instructions - they have been written to ensure you get the very best from your purchase.

Safety is Important

To ensure your safety and the safety of others, please ensure you read the Safety Instructions before you operate this product. Keep this information in a safe place for future reference.

IMPORTANT SAFEGUARDS

- 1- The principle of this device is based on the measurement of body's electronic impedance known as bio electrical term named as 'Bioelectronic Impedance Analysis (BIA), a factor related to the body fat & hydration % and in con-current with other physical data (Age, Gender, Body Height).
- 2- The evaluated readings of body fat & hydration % are only reference value. It is not recommended for use of any medical purpose. Please always consult your medical doctor for advice of any medication or diet in order to achieve normal weight level.
- 3- Pacemaker users/wearers, or other medical implants with electronic components, are not recommended to use this scale. Implant performance may be jeopardized and impaired. Pregnant and cardiac pacemaker users are prohibited to use the scale.
- 4- To prevent interference, do not use the

device, in environment near equipment with strong electro-magnetic field. Always keep a good distance from such equipments when operating the scale, or simply turn off those equipments radiating strong electromagnetic fields before operating.

- 5- This scale requires barefeet for proper conductivity while measuring body fat.
- 6- Do not place the scale on a carpeted or uneven surface.
- 7- Excessive drinking, eating, exercising, medical conditions, taking medication, the female menstrual cycle., have an effect on the measurement results.
- 8- Do not disassemble the scale as incorrect handling may cause injury.
- 9- Body fat Scale is not suitable for use people of following conditions:
  - Children under 10, and adults over 80
  - People with symptoms of fever, edema, osteoporosis.
  - People on dialysis treatment.
  - People taking cardiovascular medicine.
  - Being Pregnant
  - Athletes who undertake more than 10 hours of intensive sport activity.
  - Those with heart pulse rate under 60 beats/minute.

SPECIFICATIONS

- Body fat, water, muscle, bone, kcal analysis
- high precision strain gauge sensors
- Capacity 180kg.
- 8 mm tempered glass platform
- Data memory for up to 12 person
- LCD size: 90x50(mm) with blue backlight
- Graduation: 100g or 0.2lb,0.1 pct body fat and body water
- KG/LB options

FEMALE							
Age	Body Fat				Body Water	Muscle Propotion	Calorie (kcal) deman/day
	Under Fat	Healthy	Over fat	Obese	Normal	Normal	Suggestion
10-19	<20	20-21	31-37	>37	>56	>34	1125-2700
20-39	<21	21-32	32-38	>38	>54	>33	
40-59	<23	23-34	34-40	>40	>53	>32	
60-79	<25	25-36	36-42	>42	>52	>30	
Atlet	12-25				>56	>38	2100-4200


	MALE						
Age	Body Fat				Body Water	Muscle Propotion	Calorie (kcal) deman/day
	Under Fat	Healthy	Over fat	Obese	Normal	Normal	Suggestion
10-19	<8	8-19	19-25	>25	>58	>38	1450-3300
20-39	<9	9-20	20-26	>26	>56	>37	
40-59	<11	11-22	22-28	>28	>55	>36	
60-79	<13	13-24	24-30	>30	>54	>34	
Athlet	12-25				>60	>42	2500-5000

**Athlete Mode:** Athletic body has greater muscle mass and tend to be more dehydrated than standard adult body, and will result in high body fat reading if measured with standard adult mode. "Athlete" is defined as a person who is involved in intense physical activity of approximately 12 hours per week and who has a resting heart rate of approximately 60 beats per minute or less.

**Note:** The scope above were gained by sampling, the data is only for reference, if you have any question about your health, please consult your doctor.

### Weight Measurement Only


- Operate the scale on a secure flat floor.
- Choose pounds (lb), kilograms (kg) by switching the convertible switch
- Note:** switch should be factory-set to kg
- Press the scale, whole screen appears

- Wait until all digits light for 2 sec, and "0.0" appears.


- Stand on the scale and stand still to measure your weight.

- Step off the scale and it will shut off in 10 seconds automatically. If you stand on the scale all time, the scale will shut off in 10 seconds automatically.

### DATA SETTING

- User Setting  
Press the SET button, "-01-" will flashes and then press the UP or DOWN button to choose user number among them Press the SET button to confirm.


- **Gender Setting**  
Default gender icon will flash. Press the UP or DOWN button to choose between male and female gender.

### - Age Setting


Age appears and default age will flash  
Press the UP or DOWN button to choose  
the age and press the SET button to confirm.

### - Height Setting

Height indicator HT and inch or cm symbol  
light up (if switch set to "lb or st" default is  
height in inches;

if switch set to "kg" default is height in cm)  
Default height flashes. Press the UP or  
DOWN button to choose the height.  
Pressing and holding either the UP or  
DOWN button will advance numbers  
quickly. Press SET the button to confirm  
the height.


### - Setting Over

If SET button is pressed. User data  
displayed for 2 seconds

### BODY ANALYSIS


- Press the Up or Down button to select your  
user number among the 12 user numbers  
and recall the programmed data.


- Data appears for 2 seconds and then  
displays "0.0"

- Step on the scale barefoot with your feet  
aligned on the electrodes on the two sides  
of platform, weight displays for 3 sec, then  
test body resigistor.


- weight, "Body fat" "body water" display 3  
seconds.


- Kcal, Body muscle, Bone displays for  
seconds then turn off.

### ERROR INFORMATION

Weight overload & error in fat measuring

Replace battery

### TECHNICAL SPECIFICATION

Weight Range Max 180 kg  
Graduation 0.1 kg  
Age Range : 10 ~ 99  
Graduation 1  
Height Range : 100 ~ 220 cm  
Graduation 1 cm  
Body Fat : 5 ~ 53 %

Graduation 0.1%  
Body water : 36 ~ 85 %  
Graduation 0.1%  
Body Muscle : 25 ~ 75 %  
Graduation 0.1%  
Calorie  
Graduation 1

**CAUTION**

Keep the scale off water.

Keep the scale with care.

Don't put it on uneven floor while weighing.

Don't step on the scale edge of the scale.

Don't step on the scale with wet feet.

Measuring barefoot while weighing body fat&water percentage.

**SHINE OVERSEAS MARKETING LIMITED**

**FLAT/RM 2803 28/F CHINA RESOURCES BLDG, 26 HARBOUR RD, WANCHAI HK, P.R. CHINA**

**Tel: +86-755-86182000 Fax: +86-755-86182500**